CountAX TPC / E

Keimindikator zur Bestimmung der Gesamtkeimzahl (TPC) und dem Enterobacteriaceae-Nachweis (E) in Flüssigkeiten und auf Oberflächen. Couche nourricière d'immersion à la disposition de nombre de germe total (TPC) et pour la recherche des Enterobacteriaceae (E) dans les liquides et sur les surfaces. Indicatore di germi per la determinazione del numero totale di germi (TPC) e per la determinazione della Enterobacteriaceae (E) nei liquidi e sulle superfici. Germ indicator for the determination of the total germ count (TPC) and for the Enterobacteriaceae screen (E) from liquides and for surfaces.

Gebrauchsanweisung Notice d'utilisation Istruzioni per l'uso Instructions for use

REF AL0401

Cont. 10 Dip-Slides / Lames immergées / Piastre a immersione / Dip slides
10 Etiketten / Etiquettes / Etichette / Labels

Axon Lab AG, CH-5405 Baden
Tel. +41 (56) 484 80 80
Fax +41 (56) 484 80 99
info@axonlab.ch · www.axonlab.ch

CountAX TPC / E

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Intended Use


The CountAX Germ Indicator is used for the determination of the total germ count (TPC) and for the detect Enterobacteriaceae screen (E) from liquids and on surfaces. The germ indicator can also be used for transporting samples.

Package contents

- 10 Dip-Slides
- Labels

Basic principle

The CountAX Germ Indicator includes 2 culture media:

The TPC agar side can be used to rapidly grow the most common bacteria and fungi.

The modified VRB agar side (crystal violet, neutral red bile agar with additional glucose) can be used to grow bacteria that belong to the Enterobacteriaceae (E) family and grow in red colonies.

Composition

TPC agar	Modified VRB agar (E)
Casein peptone	Peptone
Soy peptone	Yeast extract
Sodium chloride	Sodium chloride
Lecithin	Glucose
Polysorbato 80	Bile salts
Sodium thiosulfate	Crystal violet
L-Histidine	Neutral red
Agar	Agar
pH	7.3 pH
	7.4

The culture media also contain substances that minimize the effect of disinfectants/detergents (surface tension).

Note: When samples are being taken, surfaces treated with disinfectants should not come into direct contact with the agar surfaces. If this happens, the surfaces should be treated afterwards (contaminated agar residue).

Storage

The CountAX will be stored at 15 - 25 ° C. Temperature fluctuations and draughts must be avoided. Do not freeze. Do not use the product after the expiration date (EXP) specified on the package.

Do not use the product if it is discolored or dried out, if the agar layers are detached from the plastic medium, or if there are signs of bacterial/fungal growth. Do not touch the growing colonies as they might contain pathogenic germs.

Sampling

Any contamination of the agar surfaces (growing media) must be avoided. The agar surfaces must be fully placed on the analysis material. Once the sample has been taken, screw the tube tightly shut again.

1. Unscrew the CountAX tube and remove the culture medium without touching the agar surfaces.
2. Contact inoculation

(a) The flexible holder with the bracket enables easy handling (a). Solid surfaces can be tested by firmly pressing each side of the germ indicator against the relevant surface and holding it steadily for three to four seconds (b).

(b) Dipping process

Liquid samples are tested by dipping the germ indicator into the liquid for three or four seconds until the agar surfaces are completely covered (c). Use filter paper to wipe off any excess drops from the sample.

(c) Smearing process

Semi-solid materials or objects that are hardly accessible can be tested by smearing a sterile swab carefully for example over the frame of a limited area. If the object is dry, the swab will first be moistened with sterile water. A moistened swab can also be used to obtain samples from powders (e.g. spices) or viscous fluids. After smearing the sample area, roll the swab carefully across the germ indicator's agar surfaces from left to right and from bottom to top (d).

(d) Insert the culture medium again in the tube and screw the tube tightly shut.

3. Insert the culture medium again in the tube and screw the tube tightly shut.

4. Complete the sample identification label and affix it to the tube.

5. Incubate the germ indicator tightly sealed for 1 day at 35 - 37 ° C or 2 days at 27 - 30 ° C or up to 5 days at 22 ° C.

Transporting the inoculated CountAX tubes

The inoculated tubes can be transported at 15 - 25 ° C.

Interpretation of the results


After incubation, remove the germ indicator from its tube, calculate the germ count (number of colony-forming units, CFU) and examine the color reactions.

Bacteria belonging to Enterobacteriaceae grow on the modified VRB agar (E) as red colonies. Glucose also enables some other gram-negative bacteria, e.g. Pseudomonas species to grow as red colonies. The growth of gram-positive organisms is inhibited on the modified VRB agar. The following limit values may be considered as a basis for evaluating the overall contamination level in food:

Inoculation par contact	
Clean	0 CFU / side
Contaminated	1-10 CFU / side
Highly contaminated	> 10 CFU / side

Good	< 10 ⁵ CFU/g
Permissible	10 ⁵ -10 ⁷ CFU/g
Poor	> 10 ⁷ CFU/g

The limit values for cooked foods are lower than for raw foods. The presence of Enterobacteriaceae in cooked foods always indicates wrong handling or lack of hygiene.

It is difficult to specify limit values for the total germ count on work and other surfaces since the criteria vary according to usage. However, a surface may be regarded as insufficiently cleaned if it shows 5 colonies/cm² after cleaning. This figure corresponds to about 45 colonies per slide on the TPC agar side.

The following limit values may be considered as basis for evaluating the level of contamination:

Inoculation par contact	
Clean	0 CFU / side
Contaminated	1-10 CFU / side
Highly contaminated	> 10 CFU / side

Limitation of method

When using the CountAX as a contact slide, the system is equally sensitive as the contact plate procedure.

The dipping and smearing processes have a detection limit of 1'000 CFU/ml. The permissible TPC for drinking water is too low to be detected using CountAX slides.

Only results should be compared that were obtained using the same inoculation systems.

Quality control

Growth characteristics of the medium can be tested with the following strains:

	TPC	E
Staphylococcus aureus ATCC 25923	+	-
Enterococcus faecalis ATCC 29212	+	-
Escherichia coli ATCC 25922	+	+
Pseudomonas aeruginosa ATCC 27853	+	+
Salmonella thyphimurium ATCC 14028	+	+

(+) Growth (-) Partial or complete inhibition

Waste disposal

All laboratories are responsible for proper disposal of waste in accordance with the legal requirements.

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