UrinAX CL/MC/CETRI

Transportmedium, Keimzahlbestimmung und selektive Anzucht von Bakterien im Urin
Milieu de transport, dosage des germes et développement sélectif de bactéries dans l'urine
Mezzo di trasporto, conta batterica e moltiplicazione selettiva dei batteri nell'urina
Transport medium, bacterial count and selective propagation of bacteria in urine

Gebrauchsanweisung Notice d'utilisation Istruzioni per l'uso Instructions for use

REF AL0205
BD257498 **Cont.**

- 10 Dip-Slides / Lames immergées / Piastra a immersione / Dip slides
- 10 Patientenetiketten / Etiquettes pour patient / Etichette paziente / Patient labels

In vitro diagnosticum
Usage in vitro
Per uso diagnostico in vitro
For in vitro diagnostic use

Axon Lab AG, CH-5405 Baden
Tel. +41 (56) 484 80 80
Fax +41 (56) 484 80 99
info@axonlab.ch · www.axonlab.ch

Schweiz
Axon Lab AG
Täferstrasse 15
CH-5405 Baden
Tel. +41 (0)56 484 80 80
Fax +41 (0)56 484 80 99
info@axonlab.ch
www.axonlab.ch

Deutschland
Axon Lab AG
Heinrich-Otto-Strasse 1
D-73262 Reichenbach/Stuttgart
Tel. +49 (0)7153 9226-0
Fax +49 (0)7153 9226-12
info@axonlab.de
www.axonlab.de

Suisse
Axon Lab SA
En Budron E9
CH-1052 Le Mont-sur-Lausanne
Tél. +41 (0)21 651 71 50
Fax +41 (0)21 651 71 59
info_f@axonlab.ch
www.axonlab.ch

Österreich
Axon Lab AG
Gewerbepark 1
A-6404 Polling
Tel. +43 (0)523 88 77 66
Fax +43 (0)523 88 77 66-12
info@axonlab.at
www.axonlab.at

UrinAX CL/MC/CETRI

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10

Field of application


Transport medium, bacterial count and selective propagation of bacteria in urine.

Package contents

- 10 UrinAX dip inoculum culture media
- 10 patient labels

Principle

UrinAX consists of a sealed container with 3 culture media:

The CLED culture medium is used to determine the total bacterial count in the urine. The MacConkey Agar is used for the selective propagation of Enterobacteriaceae. Cetrinide agar (CETRI) is used for the selective propagation of *Pseudomonas aeruginosa*.

Composition

Basic composition in g/l

CLED Agar	MacConkey Agar
Peptone 8.0	Peptone 20
Meat extract 3.0	Lactose 10
Lactose 10	Bile salts 1.5
L-cystein 0.128	Neutral red 0.03
Bromthymol blue 0.02	pH 7.0
pH 7.5	
CETRI Agar	
Pancreas hydrolysate from gelatin 16	Cetyltrimethyl ammonium bromide 0.3
Casein hydrolysate 10	pH 7.2

8

Storage

UrinAX should be stored at 15–25 °C. Temperature fluctuations and drafts should be avoided, and incubation should occur at 35–37 °C within 24 hours. Do not freeze. Do not use the product after the expiry date specified on the package.

Collection of urinary specimens

The urinary specimen must be inoculated into the UrinAX culture medium immediately after collection. Test results can be affected if the patient has been treated with antibiotics. In this case, the test must be performed no sooner than 48 hours after taking the last dose of medication.

Performing the test

1. Unscrew the UrinAX tube and remove the culture medium without touching the cultures.
2. Dip the culture medium into the urine mid-stream until the agar surfaces are completely covered. If the urine quantity is insufficient, pour the urine onto the agar surfaces. Allow surplus urine to run off the culture medium. Dab off the last drops of urine with filter paper.
3. Return the culture medium to the tube and close the lid.
4. Complete the patient labels and affix them to the tube.
5. Incubate the tube in an upright position in an incubator at 37 °C for 16–24 hours.

In the case of a negative result, the cultures can be incubated for a further 24 hours to detect slow-growing bacteria.


Transportation of the inoculated UrinAX

The inoculated culture medium can be transported at 15–25 °C. Incubation should occur within 24 hours.

Interpretation of results

After incubation, assess the bacteria growth on CLED Agar. Compare the bacterial count on the CLED Agar side with the following readout:

Bacterial count on CLED Agar/ml

Notes: Bacterial counts < 10⁴ bacteria/ml are usually not considered pathological. Bacterial counts between 10⁴ and 10⁵ are an uncertain result and should be repeated. Bacterial counts of > 10⁵ bacteria/ml indicate a urinary tract infection.

A mixture of different bacteria strains on the UrinAX is very probably due to contamination of the urinary specimen.

MacConkey Agar:

Assess the presence of colonies. With a positive lactose reaction, the colonies are red. With a negative lactose reaction, the colonies are colorless. Carry out biochemical or immunological identification of the bacteria.

CETRI Agar:

Assess the presence of bacteria growth. The growth of *Pseudomonas aeruginosa* germs is indicated by yellowish to blue-green colonies.

Quality control

Quality controls are carried out during the manufacture of every UrinAX lot. Growth characteristics of the medium can be tested with the following strains:

- E. coli ATCC 25922
- Staphylococcus aureus ATCC 25923
- Proteus mirabilis ATCC 12453
- Pseudomonas aeruginosa ATCC 27853

Interpret the results after a 16 to 24-hour incubation period as follows:

E. coli ATCC 25922

Yellow growth on CLED medium and pink colonies on MacConkey.

Staph. aureus ATCC 25923

Only grows on CLED medium: yellow colonies.

Proteus mirabilis ATCC 12453

Blue colouration on the CLED medium and colorless colonies on the MacConkey Agar.

Pseudomonas aeruginosa ATCC 27853

Yellowish to blue-green colonies on the 3 media CLED, MacConkey and CETRI.

Disposal

All laboratories are responsible for proper disposal of waste in accordance with the legal requirements.

References

1. Brown, V.I. et al.: J. Clin. Pathol. 18,752-756 (1965)
2. King, E.O. et al.: J. Lab. Clin. Med. 44, 301-307 (1954)
3. Handbook of Media for Clinical Microbiology 2nd Ed. 2006; Atlas, R.M. et al. CRC Taylor & Francis

9