# UrinAX CL/MC/TS

Transportmedium, Keimzahlbestimmung und selektive Anzucht von Bakterien im Urin  
Milieu de transport, dosage des germes et développement sélectif de bactéries dans l'urine  
Mezzo di trasporto, conta batterica e moltiplicazione selettiva dei batteri nell'urina  
Transport medium, bacterial count and selective propagation of bacteria in urine

## Gebrauchsanweisung Notice d'utilisation Istruzioni per l'uso Instructions for use

REF AL0201  
BD257508 Cont.

- 10 Dip-Slides / Lames immergées / Piastra a immersione / Dip slides
- 10 Patientenetiketten / Etiquettes pour patient / Etichette paziente / Patient labels

In vitro diagnosticum  
Usage in vitro  
Per uso diagnostico in vitro  
For in vitro diagnostic use

**Axon Lab AG**, CH-5405 Baden  
Tel. +41 (56) 484 80 80  
Fax +41 (56) 484 80 99  
info@axonlab.ch · www.axonlab.ch

**Schweiz**  
Axon Lab AG  
Täferstrasse 15  
CH-5405 Baden  
Tel. +41 (0)56 484 80 80  
Fax +41 (0)56 484 80 99  
info@axonlab.ch  
www.axonlab.ch

**Deutschland**  
Axon Lab AG  
Heinrich-Otto-Strasse 1  
D-73262 Reichenbach/Stuttgart  
Tel. +49 (0)7153 9226-0  
Fax +49 (0)7153 9226-12  
info@axonlab.de  
www.axonlab.de

**Suisse**  
Axon Lab SA  
En Budron E9  
CH-1052 Le Mont-sur-Lausanne  
Tél. +41 (0)21 651 71 50  
Fax +41 (0)21 651 71 59  
info\_f@axonlab.ch  
www.axonlab.ch

**Österreich**  
Axon Lab AG  
Gewerbepark 1  
A-6404 Polling  
Tel. +43 (0)523 88 77 66  
Fax +43 (0)523 88 77 66-12  
info@axonlab.at  
www.axonlab.at

# UrinAX CL/MC/TS

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## Literatur

NCCLS Publication: Quality Assurance Standards for Commercially Prepared Microbiological Culture Media; Volume 24: 14, 2004.

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
6

## Field of application


Transport medium, bacterial count and selective propagation of bacteria in urine.

## Package contents

- 10 UrinAX dip inoculum culture media
- 10 patient labels

## Principle

UrinAX consists of a sealed container with 3 culture media:

The CLED culture medium is used to determine the total bacterial count in the urine. The MacConkey Agar is used for the selective propagation of Enterobacteriaceae. The agar with TS is used for testing sensitivity to this substance.

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

7

## Interpretation of results

After incubation, assess the bacteria growth on CLED Agar. Compare the bacterial count on the CLED Agar side with the following readout:

Bacterial count on CLED Agar/ml

Notes: Bacterial counts  $< 10^4$  bacteria/ml are usually not considered pathological. Bacterial counts between  $10^4$  and  $10^5$  are an uncertain result and should be repeated. Bacterial counts of  $> 10^5$  bacteria/ml indicate a urinary tract infection.

A mixture of different bacteria strains on the UrinAX is very probably due to contamination of the urinary specimen.

MacConkey Agar:

Assess the presence of colonies. With a positive lactose reaction, the colonies are red. With a negative lactose reaction, the colonies are colorless. Carry out biochemical or immunological identification of the bacteria.

TS Agar:

Assess the presence of bacteria growth. Bacteria growth on this agar indicates a resistance to trimethoprim-sulfamethoxazole (TS).

## Quality control

Quality controls are carried out during the manufacture of every UrinAX lot. Growth characteristics of the medium can be tested with the following strains:

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